







Product Type: **Bodyjets**
 Product Code: **E044246**
 Product Name: **Wave Bodyjet**
 Description: Bodyjet con getto basculante
Bodyjet with tilting wave jet
 Dimensions: 250x80 mm
 Material: Acciaio | Stainless Steel





Note: connessioni | *connections* **G1/2" / NPT**
 Notes: **completamente ispezionabile** | *completely inspectionable*


Functions:  1 Way  Wave

Finishes: **Stainless Steel**  **316 AISI** | **316 PVD Brushed**
 High Brass **316 HB**  Copper **316 CO**  Gun Metal **316 GM**  Absolute Black **316 AB**

DIMENSIONS



FLOWRATE



PRODUCT CARE

To avoid damages to the purchased equipment, please observe the following recommendations for use:

- Clean the plumbing before connecting the device.
- Periodically clean the water blade to avoid any limestone deposits that could compromise the good functioning.

Advice: The installation of a filter and water-softener system improves and extends the life of the devices and of the plant itself.

CLEANING OF METAL PARTS:

For an easy cleaning of the metal parts of the product, simply use a soft cloth moistened in water and soap or an anti-limestone detergent.

To avoid spots caused by calcium in water, we advise drying with a cloth after each use.

Advice: Please use only natural soap-based detergents avoiding those with alcohol, hydrochloric acid or phosphoric acid, bleach, chlorine etc. and those containing abrasive powders.

CERTIFICATIONS



PACKING



1 pc